TECHNICAL SHEET



Product presentation

Professional bait gun for gel application for effective and economical treatments against cockroaches and ants.

Strengths

- Metal bait gun, solid and qualitative.
- Complete kit with holster, precision needles and user guide.

- Allows precise deposition of the gel thanks to the negative pressure system, the 5 dosing positions and the different precision needles.

Application and dose of use

Suitable for standard 30-35g gel cartridges.

ADDICT[®] BAIT GUN operating by negative pressure, during the first use it is necessary to prime the gun several times, until the gel comes out.

In order to have optimal dosing performance, be sure to test the viscosity and density of the gel used before setting the bait gun.

By adjusting the wheel, 5 different doses can be selected to fit the right dose of application:

Position	Dosage
1	0,05 g
2	0,1 g
3	0,15 g
4	0,17 g
5	0,19 g

When applying ADDICT[®] gel, we recommend positioning the wheel in position 1 and using a precision needle marked "15".



Reference & Packaging

M0119 - ADDICT Bait gun / Kit

The ADDICT kit includes:

- a metal applicator gun
- a storage holster
- a plastic tube containing 4 precision needles and replacement seals
- an user guid

Special precautions

- Keep out of reach of children.
- Never place the gun's metal nozzle in electrical outlets.
- Follow the manufacturer's instructions for the precautions when using the gel cartridge.



