# GUN KIT GEL APPLICATOR

NEGATIVE PRESSURE APPLICATOR GUN



### **Packaging**

Ref. M0138: Gel applicator gun KIT

- The kit includes:
- 1 applicator gun
- 1 storage case
- 1 plastic tube containing 5 plastic tips and 2 hoods
- 1 maintenance kit containing 1 O-ring, 1 valve inlet check valve and 1 outlet check valve
- 1 instruction manual

## Safety instructions

Keep out of reach of children.
Follow the manufacturer's instructions regarding precautions for use of the gel cartridge.

#### **Presentation**

Professional gel application gun for effective and economical treatments against cockroaches and ants.

#### Strengths

- Complete kit with case, precision needles and instructions.
- Allows precise gel deposition thanks to the negative pressure system, 5 dosing positions and different precision needles.

#### Instructions for use

Suitable for standard 30 – 35g gel syringes.

The **gel applicator gun** works by negative pressure, when using for the first time it is necessary to prime the gun several times, until the gel comes out.

For optimal dispensing performance, be sure to test the viscosity and density of the gel being used before adjusting the gun setting.

By adjusting the dial, 5 different doses can be selected to suit the correct application dose:

e
g
g
g
g
g

To apply **ADDICT®**, we recommend positioning the dial in position 1 and using a precision needle marked "15".