

# EMEROD

## PRESERVATIVE LIQUID

FOR MULTI-CATCH RODENT TRAPS

### Presentation

**EMEROD PRESERVATIVE LIQUID** is a concentrated liquid deodorizer designed for use in multi-catch rodent traps.

### Strengths

- Compatible with all rodent capture devices.
- Creates an odor barrier.
- Bacteriostatic properties.
- Preserves captured rodents for up to 3 months.
- Free from food allergens.

### Instructions for use


Pour 1L of **EMEROD PRESERVATIVE LIQUID** into the recovery tank of the trap and add 2L of water. Adjust the quantities based on the size of the trap so that the solution forms a 2 cm layer in the recovery tank.

When used with the CAPTURA multi-catch rat trap (R3246):  
Mix 5 L of **EMEROD PRESERVATIVE LIQUID** + 10 L of water in the collection tray.

Once captured, rodents can be preserved for up to 3 months.  
It is recommended to replace the liquid every 6 months in case of captures. If no rodents are captured, the liquid can be kept in the tank until its expiration date.



Made in France



### Packaging

5L can

### Precautions for use

Wear gloves and wash hands after each use.  
Do not dispose of residues in sewers or waterways.  
Do not reuse the packaging.  
Keep out of reach of children and domestic animals.  
For first aid instructions, contact the nearest poison control center. In case of medical consultation, keep the container or label available.

### Disposal

The packaging and product must be disposed of as hazardous waste under the sole responsibility of the waste holder.